



- Electrical Safety Tests Overview
- Medical Device Testing Requirements per IEC 60601-1 3rd edition
- Medical Device Testing Application Overview

Webinar Notes

Please use the Q & A utility to ask us any questions concerning the material being presented.

You can find a recording of this webinar and presentation on our Webinar Archive page, www.arisafety.com/webinars/archived-webinars/

Please contact Brittany Socha – on the chat line or email Brittany.socha@ikonixusa.com if you have any connection issues.





Bishan Patel

Presenter
Applications Engineer



Anthony Arroyo

Presenter
Applications Engineer



Syed Abidi

Panelist
Applications Leader




Brittany Socha

Organizer
Communications Leader



Learning Objectives



Common Medical Electrical Safety Tests



AC Hipot



DC Hipot



Ground Bond



Leakage
Current

What is Good Ground?

National Electrical Code

- NFPA 70
- Protect people and property from electrical hazards

NEC 250-45

- "Any exposed non-current carrying metal parts of cord & plug connected equipment which may become energized shall be grounded."

NEC 250-51

- Permanent and Continuous
- Capacity to conduct fault current
- Low impedance to limit voltage to ground.

NFPA (National Fire Protection Agency) stipulates the NEC which is adopted in all 50 U.S. States. The NEC gives requirements for grounding products and installations. The NEC defines a "good" ground.

Ground Bond Test

40A



GB test verifies the integrity of the ground connection between exposed metal and ground wire of the power cord.



High current is injected into the ground pin of the product's power cord which flows through the chassis and is returned on an accessible grounded part.




Determines if the safety ground wire is capable of handling excessive current flow in case a fault occurs and the product's insulation fails.

Ground Bond Test

40A



Example Circuit for Ground Bond Test



Ground Bond Test


40A



Ground Bond test is commonly considered a Type test and is generally performed *before* the Hipot test.

The results of this test are displayed in Ohms (Ω).

The ground conductor of a product must have a low enough impedance to handle any fault current.



Ground Bond Test

40A



Test parameters for Ground Bond test vary from standard to standard.



Manufacturers must consult the safety standard which they are trying to comply with before setting test parameters.

IEC/UL 60601-1 3rd Edition

8.6.4 Impedance and Current Carrying Capability (AC Ground Bond Test)

REQUIREMENT	PASS CRITERIA
Current = 25 A OR 1.5 * highest rated current (whichever is greater $\pm 10\%$) passed through protective earthing circuit. Frequency = 50 or 60 Hz, no load voltage ≤ 6 V	Impedance protective earthing circuit on the DUT ≤ 100 m Ω For DUTs with non-detachable supply cord, impedance for DUT ≤ 200 m Ω

Dielectric Withstand Test



Also referred to as Hipot test, it is used to determine whether the insulation of a product is able to withstand an over-voltage condition for a period of time without breaking down.



It is a deliberate application of high voltage potential between the mains input and any exposed dead-metal.



The resulting leakage current (due to the application of high voltage) is measured to determine whether a product's insulation is able to withstand the high voltage without breaking down.



This test verifies that the insulation of a product is capable of protecting the user from any leakage currents as a result of an electrical fault within the product.


Dielectric Withstand Test



Can be a Type test or a Routine test.



Used to detect possible defects in the design of a product and workmanship defects such as inadequate creepage and clearance distances.



Dielectric Withstand Test



Test is performed on both Class I and Class II products.



Test can be performed in both AC and DC mode based on the safety standard, however AC hipot test is more stressful than the DC hipot test.



Test voltage and trip settings must be specified by the manufacturer in accordance with the safety standard.

Alternatively, a d.c. test voltage equal to the peak value of the a.c. test voltage may be used.

The test voltage, for the WORKING VOLTAGE to which the insulation is subjected is greater than or equal to the value specified in Table 6.

- b) During the test, **breakdown** constitutes a failure. Insulation breakdown is considered to have occurred when the current which flows as a result of the application of the test voltage rapidly increases in an uncontrolled manner, that is, the insulation does not restrict the flow of the current. Corona discharge or a single momentary flashover is not regarded as insulation breakdown.

Means of Protection

MOP is divided into two categories

MOOP- Means Of Operator Protection

MOOP is applied to situations where the equipment will not come into contact with the patient at all. For example, IVD (in-vitro diagnostic) devices which are used only in medical laboratories, such as centrifuges, are often classified as requiring MOOP only

MOPP- Means Of Patient Protection

Any equipment that is used in a medical environment, in an area which means patients could come into contact with it, requires at least one means of patient protection (MOPP). This includes everything from hospital beds and lighting to ultrasound scanners and MRI machines, to defibrillators and dialysis machines. The highest level of protection specified in the standard requires at least two MOPPs.

Dielectric strength test requirements differ between these two categories with MOOP being more liberal.

Voltage requirements are listed in table 6 of IEC 60601-1 "AC test voltages in V rms"

Test Voltages For Solid Insulation MOP

[illegible]

It's up to the manufacturer to decide whether the equipment only requires MOP for the operator (MOOP), or whether the more stringent patient levels are required (MOPP). If the manufacturer decides MOOP is enough, it will have to back up its reasoning with a risk assessment as per ISO 14971, which examines how likely it is that a patient will come into contact with the equipment.

3rd EDITION REQUIREMENTS BY CLASSIFICATION

Classification	Isolation	Creepage	Insulation
One MOOP	1500 V ac	2.5mm	Basic
Two MOOP	3000 V ac	5mm	Double
One MOPP	1500 V ac	4mm	Basic
Two MOPP	4000 V ac	8mm	Double

Leakage Current Test



Line Current Test are performed on electrical products to measure the leakage current which could flow through a person while the product is operating.



A measuring device (MD) is used to simulate the impedance of the human body under different conditions depending upon the application of the product.



Test is run under both normal and single fault conditions and reversed polarity on the input line power at 110% of rated input.




Most commonly performed on medical equipment.

Leakage Current Test




Typical circuit for Leakage Current test



Leakage Current Test



The MD can vary from standard to standard.




$$\begin{aligned} R_1 &= 10\text{k}\Omega \pm 5\%^{1)} \\ R_2 &= 1\text{k}\Omega \pm 1\%^{1)} \\ C_1 &= 0.015\mu\text{F} \pm 5\%^{1)} \end{aligned}$$

¹⁾Non-inductive components

²⁾Impedance » measuring impedance Z


Measuring device for EN
60601-1 Medical Electrical
Equipment.



Equivalent to the above in subsequent figures.

Medical Device Testing Requirements

- IEC 60601-1 3rd Edition – Medical Electrical Equipment
 - Collateral Standards (60601-1-X)
 - Particular Standards (60601-2-X)
- Safety Tests Requirements
- Type and Routine Testing
 - Impedance and Current Carrying Capabilities (AC Ground Bond Test)
 - Dielectric Strength (Hipot Test)
 - Leakage Current and Patient Auxiliary Current (Leakage Current Test)



Ground Bond Test

40A



8.6.4 Impedance and Current Carrying Capability (AC Ground Bond Test)

REQUIREMENT	PASS CRITERIA
Current = 25 A OR 1.5 * highest rated current (whichever is greater $\pm 10\%$) passed through protective earthing circuit. Frequency = 50 or 60 Hz, no load voltage ≤ 6 V	Impedance protective earthing circuit on the DUT ≤ 100 m Ω For DUTs with non-detachable supply cord, impedance for DUT ≤ 200 m Ω

- 25 A of current passed through the ground circuit
- Resistance of the ground circuit to be ≤ 100 m Ω or ≤ 200 m Ω
- Frequency of 50 or 60 Hz

Dielectric Withstand Test



8.8.3 Dielectric Strength (Hipot Test)

REQUIREMENT	PASS CRITERIA
Test voltage - Refer to tables 6 and 7 Test time = 10 sec ramp up, 60 sec dwell & 10 sec ramp down Tested at 50Hz, 60Hz or DC equivalent ($1.414 \times$ AC test voltage)	No dielectric breakdown If DUT enclosure is non-conductive, use metal foil as conductive medium for return point.

- Test voltage depends on two factors
 - Peak Working Voltage
 - Means of Protection, MOP (Operator) or MOPP (Patient)
- Dielectric breakdown is considered a failure
- For example, 240 V mains input requires 1500 V AC Test voltage

Leakage Current Test



8.7 Leakage Current and Patient Auxiliary Current* (Leakage Current Test)

REQUIREMENT	PASS CRITERIA
Tested in Normal Conditions (NC) & Single Fault Conditions (SFC) Tested with supply at 110% highest rated mains voltage Tested at highest rated supply frequency	Earth Leakage Current ≤ 5 mA (NC) or 10 mA (SFC) Touch Current ≤ 100 μ A (NC) or 500 μ A (SFC) Patient Leakage Current - Refer to Tables 3 and 4

- Normal and Single Fault Conditions
- Mains supply at 110% voltage
- Placement of the Measuring Device (MD)

Leakage Current Test



8.7.4 Leakage Current Measurements

Earth Leakage Current: Figure 13, Section 8.7.4.5

Touch Current (Enclosure Leakage): Figure 14, Section 8.7.4.6

Patient Leakage General: Section 8.7.4.7

Patient Auxiliary (Patient Lead to Lead): Fig. 19, Section 8.7.4.8

Patient Lead to Earth: Fig. 15, Section 8.7.4.7a

Mains on Applied Part (for F-Type Patient Leads): Fig. 16, Section 8.7.4.7b

Mains on Signal I/O Ports: Fig. 17, Section 8.7.4.7c

Mains on Non-Protectively Earthed Chassis Point: Fig. 18, Section 8.7.4.7d

- Example of different type of leakage current test called out in the standard
- Required test vary depending on the type of medical device and inclusion of any patient applied parts

Leakage Current Test



Table IV


***Allowable values of continuous LEAKAGE and PATIENT AUXILIARY CURRENTS, in milliamperes**

Current	Type B		Type BF		Type CF	
	N.C.	S.F.C.	N.C.	S.F.C.	N.C.	S.F.C.
EARTH LEAKAGE CURRENT general	0,5	1 ¹⁾	0,5	1 ¹⁾	0,5	1 ¹⁾
EARTH LEAKAGE CURRENT for EQUIPMENT according to notes 2) and 4)	2,5	5 ¹⁾	2,5	5 ¹⁾	2,5	5 ¹⁾
EARTH LEAKAGE CURRENT for EQUIPMENT according to note 3)	5	10 ¹⁾	5	10 ¹⁾	5	10 ¹⁾
ENCLOSURE LEAKAGE CURRENT	0,1	0,5	0,1	0,5	0,1	0,5
PATIENT LEAKAGE CURRENT according to note 5)	d.c.	0,01	0,01	0,05	0,01	0,05
	a.c.	0,1	0,1	0,5	0,01	0,05
PATIENT LEAKAGE CURRENT (MAINS VOLTAGE on the SIGNAL INPUT PART OF SIGNAL OUTPUT PART)	—	5	—	—	—	—
PATIENT LEAKAGE CURRENT (MAINS VOLTAGE on the APPLIED PART)	—	—	—	5	—	0,05
PATIENT AUXILIARY CURRENT according to note 5)	d.c.	0,01	0,01	0,05	0,01	0,05
	a.c.	0,1	0,1	0,5	0,01	0,05

Leakage Current Test



The MD can vary from standard to standard.




$$\begin{aligned} R_1 &= 10\text{k}\Omega \pm 5\%^{1)} \\ R_2 &= 1\text{k}\Omega \pm 1\%^{1)} \\ C_1 &= 0.015\mu\text{F} \pm 5\%^{1)} \end{aligned}$$

¹⁾ Non-inductive components

²⁾ Impedance » measuring impedance Z

Measuring device for EN
60601-1 Medical Electrical
Equipment.



Equivalent to the above in subsequent figures.

Earth Leakage




Figure 3 Earth Leakage Test

Enclosure Leakage





Figure 4 Enclosure Leakage Test

Patient Leakage


60601-1 © IEC:2005

– 179 –




IEC 2398/05


Additional Tests




Leakage current



Mains Voltage on
Applied Part
Figure 16




External voltage
on SIO
Figure 17



External voltage on
metal accessible part
that is not PE
Figure 18


MOAP



Mains on SIP/SOP

60601-1 © IEC:2005

- 183 -



Medical Device Testing Application Overview



This Application note provides an example of how to test a medical device per IEC/UL 60601-1 3rd Edition using the Associated Research MedTEST System.

Table 1 lists all the instruments shown in the test setup used as an example in this document.

ction.

Table 1

APPLICATIONS CONSULTING

A DIVISION OF  IKONIX USA

We can help you ensure your company is meeting industry standard testing requirements while minimizing risks for your employees in the lab and on the production line.



Operator Education

Employ best practices to comply with agency required regulations and train operators on proper testing protocols.



Safe Workstation

Ensure work areas are safe to avoid fines, litigation, and most importantly, prevent injury or fatality.



System Validation

Create IQ/OQ documentation to validate your test system and ensure that you pass your next regulatory audit.

Ready to Improve Your Production Line?

Visit ikonixusa.com/consulting Contact applications@ikonixusa.com

Contact Us

You can find a recording of this webinar and presentation on our Webinar Archive page, arisafety.com/webinars/archived-webinars/

Check out our website for more information arisafety.com

For any questions about our webinars send an email to Amanda at brittany.socha@ikonixusa.com

